

# PCS<sup>®</sup>2

## Your choice for plasma collection



### Easy to use

- Self-loading pumps
- Auto-priming
- Comprehensive error and diagnostic messages
- Advanced optical sensors
- Easy-to-read display
- Minimal alarms
- Average procedure time  
37 min (for 600 ml)

### Multiple Safety Features

- Redundant air detectors
- Donor line pressure monitor
- Low AC ratios
- Fluid sensor
- Enclosed centrifuge well

### Multiple Protocols

- Platelet Poor Plasma
- Platelet Poor Plasma with optional saline replacement
- Source Leukocytes/Plasma

### Compact Size, Lightweight

- Easy to move
- Lightweight
- Ideal for centers with restricted space and for mobile collections





*"The PCS®2 is easy to use, set up, and monitor, which makes things easy for training."*

Lead LVN, Fresno, CA

## Specifications


PCS2 alone (closed)	height	17.375 in. (44 cm)
	width	21.5 in. (55 cm)
	depth	13 in. (33 cm)
PCS2 with cart (open)	height	55 in. (140 cm)
	width	22.25 in. (57 cm)
	depth	26.5 in. (67 cm)
	pole height	26.5 in. (from top deck)
Weight of PCS2	56 lb (25 kg)	
Weight of cart	29.2 lb (13.2 kg)	
Electrical voltage	100/220 VAC	
Operating frequency	47–63 Hz	
Current	5 amps domestic/F2.5 amps international	
Power cord length	15'9"	

For detailed information describing intended use, warnings, precautions and contraindications, refer to the instructions provided with each device or contact Haemonetics Corporation.

## Ordering Information

Description	List No.	Units/Case
Harness	620	100
Blow-molded Bowl	625B	25
Collection Bag	690	100
Collection Bag with Saline Adapter	692	100
Saline Adapter	695	100
Effluent Line Adapter	697	100

## PCS2 Disposable Setup with Saline Return



# HAEMONETICS®

HAEMONETICS CORPORATION  
 400 Wood Road, Braintree, Massachusetts 02184 USA  
 (800) 225-5297 / (781) 848-7100  
[www.Haemonetics.com](http://www.Haemonetics.com)